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<b>P1 CAL</b>	[Calibration]		
P1.1 iCAL			[internal calibration]
P1.2 ECAL			[external calibration]
P1.3 tCAL			[temperature calibration]
P1.4 ACAL		both/nonE/tenno/tinne	[automatic calibration]
P1.5 CALt		1 h = 12h	[time of automatic calibration]
P1.6 CALr		YES/no	[report from calibration]
<b>P2 GLP</b>	[Good Laboratory Practice]		
P2.1 USr		-	[name of user]
P2.2 PrJ		-	[name of project]
P2.3 PtiN		YES/no	[printout of measurement time]
P2.4 PdAt		YES/no	[printout of measurement date]
P2.5 PUSr		YES/no	[printout of user name]
P2.6 PPrJ		YES/no	[printout of project name]
P2.7 Pld		YES/no	[printout of factory number of balance]
P2.8 PFr		YES/no	[printout of frames]
<b>P3 rEAd</b>	[Main user parameters]		
P3.1 AuE		Stand/Slow/Fast	[filtering level]
P3.2 ConF		FASt_rEL/Fast/rEL	[confirmation of result]
P3.3 Auto		On/OFF	[autozero]
P3.4 Ldi0		ALuuAYS/never/uu_StAb	[light off the last digit]

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
<b>P4 Print</b>	[Data transmission – RS 232]		
P4.1 bAud		2400/4800/9600/19200	[speed of transmission]
P4.2 CntA		YES/no	[continuous printout in basic unit]
P4.3 Cntb		YES/no	[continuous printout in currently used unit]
P4.4 rEPL		YES/no	[manual or automatic operating]
P4.5 PStb		YES/no	[result printout: stable or unstable]
P4.6 Lo		000.0000	[min mass for automatic operating]
<b>P5 Unit</b>	[Mass units]		
P5.1 StUn		g/mg/ct/oz/ozt/dwt/mom/G	[basic mass unit – e.g. gram]
P5.2 mg		YES/no	[mg – milligram]
P5.3 Ct		YES/no	[ct – carat]
P5.4 oZ		YES/no	[oz – ounce]
P5.5 oZt		YES/no	[ozt – troy ounce]
P5.6 dwt		YES/no	[dwt – pennyweight]
P5.7 t		YES/no	[t – tael]
P5.8 nno		YES/no	[mom – momme]
P5.9 Gr		YES/no	[gr – grain]
<b>P6 Func</b>	[Functions]		
P6.1 FFun			[choice of balance function]
P6.2 PcS		YES/no	[counting pieces]
P6.3 HiLo		YES/no	[check weighing]
P6.4 PrcA		YES/no	[% deviation with reference to weighed standard mass]
P6.5 ProB		YES/no	[% deviation with reference to declared standard mass]
P6.6 d_Co		YES/no	[density determination of solids]
P6.7 d_Li		YES/no	[density determination of liquids]
<b>P7 othEr</b>	[Other functions]		
P7.1 bL		On/Aut/OFF	[display backlight]
P7.2 bEEP		On/OFF	[keys sound]
P7.3 PmS			[printout of balance parameters]

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